



## Rabbit anti-Goat IgG Antibody

<b>Product Code</b>	SGJH004
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze
<b>Immunogen</b>	IgG purified from Goat serum
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Goat
<b>Applications</b>	ELISA, WB; IHC
<b>Recommended dilution</b>	1:1000-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300
	Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	salt precipitation combined with affinity chromatograph