



## Goat Anti-Rabbit IgG Antibody;HRP conjugated

**Product Code** SGJH009

**Storage** Upon receipt, store at -20°C or -80°C. Avoid repeated freeze

**Immunogen** IgG purified from Rabbit serum

**Raised In** Goat

**Species Reactivity** Rabbit

**Applications** ELISA,WB;IHC

**Recommended dilution** 1:1000-1:5000

**Form** Liquid

**Conjugate** HRP-conjugated

**Storage Buffer** Preservative: 0.03% Proclin 300

Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

**Purification Method** salt precipitation combined with affinity chromatograph