



## Goat Anti-Mouse IgG Antibody;Alkaline phosphatase conjugated

**Product Code** SGJH011

**Storage** Upon receipt, store at -20°C or -80°C. Avoid repeated freeze

**Immunogen** IgG purified from Mouse serum

**Raised In** Goat

**Species Reactivity** Mouse

**Applications** ELISA,WB;IHC

**Recommended dilution** 1:1000-1:5000

**Form** Liquid

**Conjugate** Alkaline phosphatase-conjugated

**Storage Buffer** Preservative: 0.03% Proclin 300

Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

**Purification Method** salt precipitation combined with affinity chromatograph