



Goat Anti-Rabbit IgG Antibody;Alkaline phosphatase conjugated

Product Code SGJH012

Storage Upon receipt, store at -20°C or -80°C. Avoid repeated freeze

Immunogen IgG purified from Rabbit serum

Raised In Goat

Species Reactivity Rabbit

Applications ELISA,WB;IHC

Recommended dilution 1:1000-1:5000

Form Liquid

Conjugate Alkaline phosphatase-conjugated

Storage Buffer Preservative: 0.03% Proclin 300

Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Purification Method salt precipitation combined with affinity chromatograph